


MYFLON

Meat and Food Pads

Extend the shelf life and enhance the appearance of your meat and other food products. Our customizable packaging solutions provide the perfect absorption levels for each type of food.



Presentation is the key to success!

A fresh, appetizing presentation of goods signals quality and significantly influences the purchase decision. On the way from producer to consumer, Myflon pads, which are individually customized to the particular food, ensure high-quality freshness, hygiene and appearance.

Customized solutions: less food waste - more profit.



all-round

Customized for all solutions!

Customized absorbent pads from Myflon match all specific meat requirements and fulfil the challenge to absorb all free liquid within the packaging.

- Absorption is measured in ml/m²
- Available in rolls or already cut pads
- Liquid absorption is achieved through the surface via the micro-pore technology
- Absorption of exudate through the open sides (also in tilted position)b

MeatPad

MP - our all-round performance pad

Product range examples:	
high	MP 4000 MP 2500
medium	MP 1700 MP 1500 MP SL
low	MP 801

- For low to high absorption*
- Ideal for MAP packaging
- Available in various grades and coloursb

* Test method: measured with 0,2% NaCl and 0,1 Psi



high performance

High performance under pressure.

Due to higher absorption, MeatGuard is well suited for products that release large amount of exudate. The liquid is bonded within the core and retained there even under pressure.

- Absorption is measured in ml/pad
- Available in rolls, already cut pads or festooned
- Liquid absorption is achieved through the 3D-surface-system
- Exudate is bound in the core and will not leak from the padb

MeatGuard

MG - our high performance pad

Product range examples:	
very high	MG 150-L MG 90-S MG 90-L
high	MG 75-S MG 75-L MG 40-S

- For high to very high absorption*
- Ideal for MAP packaging
- Available in three colours and two formats
- Also available as cross productionb

* Test method: measured with 0,2% NaCl and 0,1 Psi

L= 102mm x 170mm

MYFLON

Eskisehir Organize san. Bol. 26. Cad. No 9 Eskisehir , Turkey

www.mulcfilm.com felix@felixplastic.com

